


# User Manual

## 1, Parts of Trolley


### (1). Front




(2). Back



2、 Specification





Rotating Arm Length : 1100mm

Balance Arm Length : 950mm

Lamp Head Max Width : 790mm

Max Total Length : 2470mm

3、 Installation


Unpack and take out the below 3 main parts:



1. Trolley



2. Rotating Arm



3. Balance Arm + Lamp Head



(1)、 Take off the circlip which is at the connection between rotating arm and trolley top. (Tool Needed : circlip pliers)



(2)、 Insert rotating arm(Loosen damping screw)



(3)、 Put on the circlip of rotating arm




(4)、 Pass the bundle of cables through rotating arm



(5)、 Take off the circlip which is at the connection between rotating arm and balance arm



6、 Match the number symbols and connect terminals, then install balance arm



7、 Put on the circlip of balance arm



8、 Put on the screw



9、 Install the monitor



#### 4、Warning

Please notice that the two side of balance arm should be matching as below images :

